

SWAT-T™

STRETCH-WRAP-AND-TUCK TOURNIQUET



RESEARCH

2015: Different Width and Tightening System: Emergency Tourniquets on Distal Limb Segments.

Journal of Special Operations Medicine
Conclusion: SWAT-T™ vs. CAT/RMT/SOFTT-W Better Occlusion, Safer Pressures; 10.4% failure in other TQs, zero failures in SWAT-T™ group.

2014: Assessment of 32,956 SWAT-T™ applications.

TEMS Solutions: August 2014
Conclusion: The SWAT-T™ is a durable medical device with zero first-time-use failures, when properly used.

2013: Tourniquets and occlusion: the pressure of design.

Military Medicine: May 2013
Conclusion: The SWAT-T™ performed better than the CAT (Combat Application Tourniquet). The SWAT-T™ had safer pressures, and was more effective.

2012: Lighting did not affect self-application of a stretch and wrap style tourniquet.

Journal of Special Operations Medicine
Conclusion: The SWAT-T™ stretch and wrap style tourniquet can be self-applied properly even in darkness. When properly applied, it stops limb arterial flow.

2012: Stretch and wrap style tourniquet effectiveness with minimal training.

Military Medicine: November 2012
Conclusion: The SWAT-T™ can easily be properly applied and can stop arterial flow at a variety of extremity locations. Proper application is associated with cessation of arterial flow.

Innovation and Versatility in **TRAUMA CARE**

In a class by itself, the SWAT-T™ is a multi-function trauma care device. Multiple Documented Saves!



Designed by a former SOF Operator/Medic, the SWAT-T will treat a variety of injuries, minor to life-threatening. The SWAT-T™ is being carried by Military (Conventional and Special Operations-medical and non-medical personnel), EMS, Law Enforcement Officers, Contractors, and Federal Agents. Carried by many as a pressure dressing, all-purpose wrap, primary and/or back-up tourniquet.

"The ultimate trauma multi-tool." -HA US Navy, SEAL

"The SWAT-T packs light and is superior for kids." -RS, Special Forces 18D

"The KISS principle on steroids!" -PT, Operations Manager, Triple Canopy

"3 SWAT-Ts used after Boston blast worked amazingly well." -NM

"Worked great after windlass failure (high axillary wound)."

-CY - US Border Patrol

"The SWAT-T saved my K9 partner! As a tactical team member, I will never deploy without a SWAT-T." -NL

"It was soaked in blood and it held up great, still plenty of friction. Successful application, bleeding terminated." -SV



SWAT-T™ TRAINING VIDEOS@:
SWAT-T.com



GSA Available
 NSN Pending
 Patents Pending

